


Product Name: Recombinant Human FOXD4L1 Protein
Catalog #: PEH5348



Summary

Name	FX4L1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NU39
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DEEEEASQKF LEQSLQPLQ VARWGGVALP REHIEGGGPS DPSEFGTEFR APPRSAASE DARQPAKPPY SYIALITMAI LQSPHKRLTL SGICAFISGR FPYRRKFPA WQNSIRHNLS LNDCFVKIPR EPGHPGKGTY WSLDPASQDM FDNGSFLRRR KRFKRHQLTP GAHLPHPFPL PAAHAALHNP RPGLLGAPA LPQPVPAYP NTAPGRRPYA LLHPHPPRYL LLSAPAYAGA PKKAEGADLA TPGTLPLVLP SLGPQPWEEG

image



Background

Alternative Names Forkhead box protein D4-like 1 (FOXD4-like 1)

Product Name: Recombinant Human FOXD4L1 Protein
Catalog #: PEH5348



Note

For Research Use Only , Not for Diagnostic Use