


## Summary

|                              |   |
|------------------------------|---|
| <b>Name</b>                  | SDA1_HUMAN  |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE   |
| <b>Endotoxin level</b>       | Please contact with the lab for this information  |
| <b>Tag</b>                   | C-His   |
| <b>Accession #</b>           | Q9NVU7  |
| <b>Host</b>                  | E.coli  |
| <b>Species</b>               | Human   |
| <b>Molecular Mass</b>        | 38 kDa  |
| <b>Buffer</b>                | Lyophilized from Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.  |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.   |
| <b>Sequence</b>              | TKILLFAAQA SHHLVPPEII QLLMTVANN FVTDKNSGEV MTVGINAIKE<br>ITARCPLAMT EELLQDLAQY KTHKDKNVMM SARTLIHLFR TLNPQMLQKK<br>FRGKPTIASI EARVQEYGEI DAKDYIPGAE VLEVEKEENA ENEDGWEST<br>SLSEEDADG EWIDVQHSSD EEQQEISKKL NSMPMEERKA KAAAISTSRV<br>LTQEDFQKIR MAQMRKELDA APGKSQKRKY IEIDSDEEPR GELLSLRDIE<br>RLHKKPKSDK ETRLATAMAG KTDRKEFVRK KTKTNPFSST TNKEKKKQKN<br>FMMMRYSQNV |

## image



## Background

|                          |  |
|--------------------------|--|
| <b>Alternative Names</b> | Protein SDA1 homolog (Nucleolar protein 130) (SDA1 domain-containing protein 1) (hSDA) |
|--------------------------|--|

**Product Name: Recombinant Human SDAD1 Protein**  
**Catalog #: PEH5356**

---



**Note**

For Research Use Only , Not for Diagnostic Use