


Product Name: Recombinant Human IFT57 Protein
Catalog #: PEH5358



Summary

Name	IFT57_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NWB7
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGLEDGVPRS RGEPTGEVVL ERGPGAAHYM FVVMEDLVEK LKLLRYEEEF LRKSNLKAPS RHYFALPTNP GEQFYMFCTL AAWLINKAGR PFEQPQEYDD PNATISNILS ELRSFGRTAD FPPSKLKSGY GEHVCYVLDC FAEELKYIG FTWKRPIYPV EELEESVAE DDAELTLNKV DEEFVEEETD NEENFIDLNV LKAQTYHLDL NETAKQEDIL ESTTDAAEWS LEVERVLPQL KVTIRTDNND WRIHVDQMHQ HRSIESALK ETKGFLDKLH NEITRTLEKI SSREKYINNQ LENLVQEYRA AQAQLSEAKE RYQQGNGGVT ERTRLLSEVM EELEKVKQEM EEKGSSMTDG APLVKIKQSL TKLKQETVEM DIRIGIVEHT

image



Background

Alternative Names	Intraflagellar transport protein 57 homolog (Dermal papilla-derived protein
--------------------------	---

Product Name: Recombinant Human IFT57 Protein
Catalog #: PEH5358



8) (Estrogen-related receptor beta-like protein 1) (HIP1-interacting protein)
(MHS4R2)

Note

For Research Use Only , Not for Diagnostic Use