


Product Name: Recombinant Human SLTM Protein
Catalog #: PEH5359



Summary

Name	SLTM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NWH9
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEEKRSEKS EKKESKDTKK IEGKDEKNDN GASGQTSESI KKSEEKKRIS SKSPGHMVIL DQTKGDHCRP SRRGRYEKIH GRSKEKERAS LDKKRDKDYR RKEILPFEKM KEQRLREHLV RFERLRRAME LRRRREIAER ERRERERIRI IREREERERL QRERERLEIE RQKLERERME RERLERERIR IEQERRKEAE RIAREREELR RQQQQLRYEQ EKRNSLKRPR DVDHRRDDPY WSENKKLSLD

image



Background

Alternative Names	SAFB-like transcription modulator (Modulator of estrogen-induced transcription)
--------------------------	---

Note

Product Name: Recombinant Human SLTM Protein
Catalog #: PEH5359



For Research Use Only , Not for Diagnostic Use