


Product Name: Recombinant Human ZNF446 Protein
Catalog #: PEH5360



Summary

Name	ZN446_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NWS9
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEEPLGSPHP SGTVESPGEG PQDTRIEGSV QLSCSVKEEP NVDGQEVAPS SPPLAAQSPE GNHGHQEPAS TSFHPPRIQE EWGLLDRSQK ELYWDAMLEK YGTVVSLGLP PHQPEAQAS ELGMILLTGTG VCRSLRSGNE SEGPPGCPEA QPPQGP GPAA WEGLSGAATP APTVRPGTTP

image



Background

Alternative Names Zinc finger protein 446 (Zinc finger protein with KRAB and SCAN domains 20)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF446 Protein
Catalog #: PEH5360

