


Product Name: Recombinant Human FAM120C Protein
Catalog #: PEH5361



Summary

Name	F120C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NX05
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IPPQGGKLEI AGMVVGQWAG SRSSRGRGSF GMQVVSVGGP GKGHGKEQTG RGSKGHKKGN KQGSSDGVSK SLELHQGRSR SQVNGNSGAL IKEEKSDHRL PAPSQCALS RDSNECNNGNR YLPMNNREKN HLQEQKLETV

image



Background

Alternative Names Constitutive coactivator of PPAR-gamma-like protein 2 (Protein FAM120C)
(Tumor antigen BJ-HCC-21)

Note

For Research Use Only , Not for Diagnostic Use