


Product Name: Recombinant Human NSMCE4A Protein
Catalog #: PEH5363



Summary

Name	NSE4A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NXX6
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSARERREAP ERPSLEDTEP SDSDGDEMMDP ASLEAEADQG LCRQIRHQYR ALINSVQQNR EDILNAGDKL TEVLEEANTL FNEVSRAREA VLD AHFLVLA SDLGKEKAKQ LRS DLSSFD M LRYVETLLTH MGVNPLEAEE LIRDESPDF EFIVYDSWKI TGR TAENTFN KTHTFHFLG SIYGECVVPK PRVDRPRKVP VIQEERAMPA QLRRMEESHQ

image



Background

Alternative Names	Non-structural maintenance of chromosomes element 4 homolog A (NS4EA) (Non-SMC element 4 homolog A)
--------------------------	--

Note

Product Name: Recombinant Human NSMCE4A Protein
Catalog #: PEH5363



For Research Use Only , Not for Diagnostic Use