


Product Name: Recombinant Human RBAK Protein
Catalog #: PEH5370



Summary

Name	RBAK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NYW8
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDVAVDFTQE EWQQLDPDEK ITYRDVMLEN YSHLVSVGVD TTKPNVVIKL EQGEEPWIMG GEFPCQHSPE AWRVDDLIER IQENEDKHSR QAACINSKTL

image



Background

Alternative Names	RB-associated KRAB zinc finger protein (RB-associated KRAB repressor) (hRBaK) (Zinc finger protein 769)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use