


Product Name: Recombinant Human PKD2L2 Protein
Catalog #: PEH5374



Summary

Name	PK2L2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NZM6
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YSIGRRLD FE LGKMIKQSYK NVLEKFRLKK AQKDEDK KTK GSGDLAEQAR REGFDENEIQ NAEQMKKWKE RLEKKYYSME IQDDYQPVTQ EEFRELFLYA

image



Background

Alternative Names Polycystic kidney disease 2-like 2 protein (Polycystin-2L2) (Polycystin-L2)

Note

For Research Use Only , Not for Diagnostic Use