


Product Name: Recombinant Human CCDC59 Protein
Catalog #: PEH5375



Summary

Name	TAP26_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P031
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EERHRKQARK VDHPLSEQVH QPLLEEQCSI DEPLFEDQCS FDQPQPPEEQC IKTVNSFTIP KKNKKKTSNQ KAQEEYEQIQ AKRAAKKQEF ERRKQEREEA QRQYKKKKME VFKILNKKTK KGQPNLNVQM EYLLQKIQEK

image



Background

Alternative Names	Thyroid transcription factor 1-associated protein 26 (TTF-1-associated protein 26) (Coiled-coil domain-containing protein 59) (TTF-1-associated protein BR2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CCDC59 Protein
Catalog #: PEH5375

