


Product Name: Recombinant Human ZKSCAN7 Protein
Catalog #: PEH5376



Summary

Name	ZKSC7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P0L1
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GNSGDQAGAT VLRMVRPQDT VAYEDLSVDY TQKKWKSLLT SQRALQWNMM PENHHSMA SL AGENMMKGSE LTPKQEFFKG SESSNRTSGG LFGVVPGAEE TGDVCEDTFK ELEGQTSDEE GSRLNDFLE ITDEDKKKST

image



Background

Alternative Names	Zinc finger protein with KRAB and SCAN domains 7 (Zinc finger protein 167) (Zinc finger protein 448) (Zinc finger protein 64)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use