


Product Name: Recombinant Human CEP72 Protein
Catalog #: PEH5378



Summary

Name	CEP72_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P209
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TAAQDSSAMV GEDVGSLAE SKSLQSRLAE QQQQHAREMS EVTAE LHHTH KELDDL RQHL DKSLEEN SRL KSLLSMKKE VKSADTAATL NLQIAGLQTS VKRLCGEIVE LKQHLEHYDK IQELTQMLQE SHSSLVSTNE HLLQE LSQVR

image



Background

Alternative Names Centrosomal protein of 72 kDa (Cep72)

Note

For Research Use Only , Not for Diagnostic Use