


Product Name: Recombinant Human POGK Protein
Catalog #: PEH5379



Summary

Name	POGK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P215
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLKEEEEEEE IQSRELEDGP ADMQKVRICS EGGWVPALFD EVAIYFSDEE WEVLTEQQKA LYREVMRMNY ETVLSLEFPF PKPDMITRLE GEEESQNSDE WQLQGGTSAE NEESDVKPPD WPNPMNATSQ FPQPQHFDSE GLRLPRDITE LPEWSEGYPF YMAMGFPGYD LSADDIAGKF QFSRGMRRSY DAGFKLMVVE YAESTNNCQA AKQFGVLEKN

image



Background

Alternative Names Pogo transposable element with KRAB domain

Note

Product Name: Recombinant Human POGK Protein
Catalog #: PEH5379



For Research Use Only , Not for Diagnostic Use