


Product Name: Recombinant Human CPSF2 Protein
Catalog #: PEH5382



Summary

Name	CPSF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2I0
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSQVEWMSDK LMRCFEDKRN NPFQFRHLSL CHGLSDLARV PSPKVVLASQ PDLECGFSRD LFIQWCQDPK NSIILTYRTT PGTLARFLID NPSEKITEIE LRKRVKLEGK ELEEYLEKEK LKKEAAKKLE QSKEADIDSS DESDIEEDID QPSAHKTKHD LMMKGEGSRK GSFFKQAKKS YPMFPAPEER IKWDEYGEII KPEDFLVPEL QATEEEESKL ESGLTNGDEP MDQDLSDVPT KCISTTESIE IKARVTYIDY EGRSDGDSIK KIINQMKPRQ LIIVHGPPEA SQDLAECCRA FGGKDIKVYM PKLHETVDAT SETHIQVRL KDSLVSSLQF CKAKDAELAW IDGVLDMRVS KVDTGVILEE GELKDDGEDS EMQVEAPSDS SVIAQQKAMK SLFGDDEKET

image



Background

Product Name: Recombinant Human CPSF2 Protein
Catalog #: PEH5382



Alternative Names

Cleavage and polyadenylation specificity factor subunit 2 (Cleavage and polyadenylation specificity factor 100 kDa subunit) (CPSF 100 kDa subunit)

Note

For Research Use Only , Not for Diagnostic Use