


Product Name: Recombinant Human RCOR3 Protein
Catalog #: PEH5384



Summary

Name	RCOR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2K3
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEGNTEQPVQ TSKIGLGRRE YQSLQHRHHS QRSKCRPPKG MYLTQEDVVA VSCSPNAANT ILRQLDMELI SLKRQVQNAK QVNSALKQKM EGGIEEFKPP ESNQKINARW TTEEQLLAVQ GVRKYGKDFQ AIADVIGNKT VGQVKNFFVN YRRRFNLEEV LQEWAEQGT QASNGDASTL GEETKSASNV PSGKSTDEEE EAQTPQAPRT

image



Background

Alternative Names REST corepressor 3

Note

Product Name: Recombinant Human RCOR3 Protein
Catalog #: PEH5384



For Research Use Only , Not for Diagnostic Use