


Product Name: Recombinant Human RERE Protein
Catalog #: PEH5386



Summary

Name	RERE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2R6
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DCRIHFKKYG ELPIIEKPVD PPPFMFKPVK EEDDGLSGKH SMRTRRSRGS MSTLRSGRKK QPASPDGRTS PINEDIRSSG RNPSAASTS SNDSKAETVK KSAAKKVKEEA

image



Background

Alternative Names	Arginine-glutamic acid dipeptide repeats protein (Atrophin-1-like protein) (Atrophin-1-related protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use