



Product Name: Recombinant Human SUCO Protein
Catalog #: PEH5389



Summary

Name	SUCO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBS9
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEIDNETEQK SESFSSIEKP SITYETNKVN ELMDNIIKED VNMQIFTKL SETIVPPINT ATVPDNEDGE AKMNIADTAK QTLISVVDSS SLPEVKEEEQ SPEDALLRGL QRTATDFYAE LQNSTDLYA NGNLVHGSNQ KESVFMRLNN RIKALEVNMS LSGRYLEELS QRYRKQMEEM QKAFNKTIVK LQNTSRIAE QDQRQTEAIQ LLQAQLTNMT QLVSNLSATV

image



Background

Alternative Names	SUN domain-containing ossification factor (Membrane protein CH1) (Protein osteopotential homolog) (SUN-like protein 1)
--------------------------	--

Note

Product Name: Recombinant Human SUCO Protein
Catalog #: PEH5389



For Research Use Only , Not for Diagnostic Use