

Product Name: Recombinant Human HMGXB4 Protein
Catalog #: PEH5395



Summary

Name	HMGX4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UGU5
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LHTDGHSEKK KKKEEKDKER ERGEKPKKKN MSAYQVFCKE YRVTIVADHP GIDFGELSCK LAEVWKQLPE KDKLIWKQKA QYLQHKQNKA EATTVKKRAS SSEGS MKVKA

image



Background

Alternative Names HMG domain-containing protein 4 (HMG box-containing protein 4) (High mobility group protein 2-like 1) (Protein HMGBCG)

Note

For Research Use Only , Not for Diagnostic Use