


Product Name: Recombinant Human SUN2 Protein
Catalog #: PEH5396



Summary

Name	SUN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UH99
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLEGLVSRRE AALKEDFRRE TAARIQEELS ALRAEHQQDS EDLFKKIVRA SQESEARIQQ LKSEWQSMTQ ESFQESSVKE LRRLEDQLAG LQQELAALAL KQSSVAEEVG LLPQQIQAVR DDVESQFPAW

image



Background

Alternative Names	SUN domain-containing protein 2 (Protein unc-84 homolog B) (Rab5-interacting protein) (Rab5IP) (Sad1/unc-84 protein-like 2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use