


Product Name: Recombinant Human SLC24A2 Protein
Catalog #: PEH5399



Summary

Name	NCKX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UI40
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSAARDKDEP TLPAPRLQR GGSSASLHNS LMRNSIFQLM IHTLDPLAEE LGSYGKLYY DTMTEGRFR EKASILHKIA KKKCHVDENE RQNGAANHVE KIELPNSTST DVEMTPSSDA SEPVQNGNLS HNIEGAEAQT ADEEEDQPLS LAWPSETRKQ

image



Background

Alternative Names	Sodium/potassium/calcium exchanger 2 (Na(+)/K(+)/Ca(2+)-exchange protein 2) (Retinal cone Na-Ca+K exchanger) (Solute carrier family 24 member 2)
--------------------------	--

Note

Product Name: Recombinant Human SLC24A2 Protein
Catalog #: PEH5399



For Research Use Only , Not for Diagnostic Use