


**Product Name: Recombinant Human BAZ1B Protein**  
**Catalog #: PEH5400**



## Summary

|                              |                                                                                                                                                                                                                                                                                                                                                                                                         |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | BAZ1B_HUMAN                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                     |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                                                                                                                                        |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Accession #</b>           | Q9UIG0                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Molecular Mass</b>        | 40 kDa                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                      |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                          |
| <b>Sequence</b>              | <p>EQRKEYLKKK REELKKKLKE KAKERREKEM LERLEKQKRY EDQELTGKNL<br/>PAFRLVDTP E GLPNTLFGDV AMVVEFLSCY SGLLLPDAQY PITAVSLMEA<br/>LSADKGGFLY LNRVLVILLQ TLLQDEIAED YGELGMKLSE IPLTLHSVSE<br/>LVRLCLRRSD VQEESEGST DDNKDSAAFE DNEVQDEFLE KLETSEFFEL<br/>TSEELQILT ALCHRILMTY SVQDHMETRQ QMSAELWKER LAVLKEENDK<br/>KRAEKQKRKE MEAKNKENGK VENGLGKTDR KKEIVKFEPQ VDTEAEDMIS<br/>AVKSRLLAI QAKKEREIQE REMKVKLERQ</p> |

## image



## Background

|                          |                                                                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Tyrosine-protein kinase BAZ1B (EC 2.7.10.2) (Bromodomain adjacent to zinc finger domain protein 1B) (Williams syndrome transcription factor) |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|

**Product Name: Recombinant Human BAZ1B Protein**  
**Catalog #: PEH5400**



---

(Williams-Beuren syndrome chromosomal region 10 protein) (Williams-Beuren syndrome chromosomal region 9 protein) (hWALp2)

**Note**

For Research Use Only , Not for Diagnostic Use