


Product Name: Recombinant Human TRMT6 Protein
Catalog #: PEH5402



Summary

Name	TRM6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UJA5
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KVDSLLHGTF SAKMLSSEPK DSALVEESNG TLEEKQASEQ ENEDSMAEAP ESNHPEDQET METISQDPEH KGPKERGSKK DYIQEKQRRQ EEQRKRHLEA AALLSERNAD GLIVASRFHP TPLLLSLDF VAPSRPFVVY CQYKEP LLEC YTKLRERGGV INLRLSETWL

image



Background

Alternative Names	tRNA (adenine(58)-N(1))-methyltransferase non-catalytic subunit TRM6 (mRNA methyladenosine-N(1)-methyltransferase non-catalytic subunit TRM6) (tRNA(m1A58)-methyltransferase subunit TRM6) (tRNA(m1A58)MTase subunit TRM6)
--------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Product Name: Recombinant Human TRMT6 Protein
Catalog #: PEH5402



Note

For Research Use Only , Not for Diagnostic Use