


Product Name: Recombinant Human GTF3C4 Protein
Catalog #: PEH5408



Summary

| | |
|------------------------------|---|
| Name | TF3C4_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9UKN8 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 10 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | SPGMGNADDE QQEETSSKQ VVKQGLQERS KEGDVEEPTD DSLPTTG DAG GREPMEEKLL EIQGKIEAVE |

image



Background

| | |
|--------------------------|--|
| Alternative Names | General transcription factor 3C polypeptide 4 (EC 2.3.1.48) (TF3C-delta) (Transcription factor IIIC 90 kDa subunit) (TFIIIC 90 kDa subunit) (TFIIIC90) (Transcription factor IIIC subunit delta) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use