


Product Name: Recombinant Human INTS6 Protein
Catalog #: PEH5409



Summary

Name	INT6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UL03
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PVVNNHIGGK GPPAPTTQAQ PDLIKPLPLH KISSETTNSI IHDVVENHVA DQLSSDITPN AMDTEFSASS PASLLERPTN HMEALGHDHL GTNDLTVGGF LENHEEPRDK EQCAEENIPA

image



Background

Alternative Names Integrator complex subunit 6 (Int6) (DBI-1) (Protein DDX26) (Protein deleted in cancer 1) (DICE1)

Note

For Research Use Only , Not for Diagnostic Use