


Product Name: Recombinant Human ZDHC8 Protein
Catalog #: PEH5412



Summary

Name	ZDHC8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULC8
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SIQERKDREE RERLLRSQAD SLFGDSGVYD APSSYSLQQA SVLSEGPRGP ALRYGSRDDL VAGPGFGGAR NPALQTSLS LSSSVSRAPR TSSSSLQADQ ASSNAPGPRP SSGSHRSPAR QGLPSPPGTP

image



Background

Alternative Names	Palmitoyltransferase ZDHC8 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 8) (DHHC-8) (Zinc finger protein 378)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use