


Product Name: Recombinant Human SLC39A10 Protein
Catalog #: PEH5415



Summary

Name	S39AA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULF5
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YKQQRGKQKW FMKQNTTEEST IGRKLS DHKL NNTPDSDWLQ LKPLAGTDDS VVSEDRLNET ELTDLEGQQE SPPKNYLCIE EEKIIDHSHS DGLHTIHEHD LHAAAHNHHG

image



Background

Alternative Names	Zinc transporter ZIP10 (Solute carrier family 39 member 10) (Zrt- and Irt-like protein 10) (ZIP-10)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use