


Product Name: Recombinant Human KIDINS220 Protein
Catalog #: PEH5417



Summary

Name	KDIS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULH0
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RVPKSPEHSA EPIRTFIKAK EYLSDALLDK KDSSDSGVRS SESSPNHSLH NEVADDSQLE KANLIELEDD SHSGKRGIPH SLSGLQDPPII ARMSICSEDK KSPSECSLIA SSPEENWPAC QKAYNLNRTP STVTLNNSA PANRANQNFD EMEGIRETSQ VILRPSSSPN

image



Background

Alternative Names	Kinase D-interacting substrate of 220 kDa (Ankyrin repeat-rich membrane-spanning protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KIDINS220 Protein
Catalog #: PEH5417

