


Product Name: Recombinant Human GUCA1B Protein
Catalog #: PEH5425



Summary

Name	GUC1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UMX6
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DEEASQYVEG MFRAFDKNGD NTIDFLEYVA ALNLVLRGTL EHKLKWTFKI YDKDGNCGCID RLELLNIVEG IYQLKKACRR ELQTEQGQQL TPEEVVDRIF LLVDENGDGQ LSLNEFVEGA

image



Background

Alternative Names	Guanylyl cyclase-activating protein 2 (GCAP 2) (Guanylate cyclase activator 1B)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use