


Product Name: Recombinant Human PCDHB10 Protein
Catalog #: PEH5426



Summary

Name	PCDBA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UN67
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGSFVVNLAK DLGLAEGELA ARGTRVVSDD NKQYLLDSH TGNLLTNEKL DREKLCGPKE PCMLYFQILM DDPFQIYRAE LRVRDINDHA PVFQDKETVL KISENTAEGT AFRLERAQDP DGGLNGIQNY TISPNSFFHI NISGGDEGMI YPELVLDKAL DREEQGELSL TLTALDGGSP SRSGTSTVRI VVLDVNDNAP QFAQALYETQ APENSPIGFL IVKVWAEDVD SGVNAEVSYS FFDASENIRT TFQINPFSGE IFLRELLDYE LVNSYKINIY AMDGGGLSAR CRVLVEVLDT

image



Background

Alternative Names	Protocadherin beta-10 (PCDH-beta-10)
--------------------------	--------------------------------------

Product Name: Recombinant Human PCDHB10 Protein
Catalog #: PEH5426



Note

For Research Use Only , Not for Diagnostic Use