


Product Name: Recombinant Human G3BP2 Protein
Catalog #: PEH5429



Summary

Name	G3BP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UN86
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GFPPRGPRPG RGDMEQNDS NRRIIRYPDS HQLFVGNLPH DIDENELKEF FMSFGNVVEL RINTKGVGGK LPNFGFVFD DSEPVQRILI AKPIMFRGEV RLNVEEKTR AARERETRGG GDDRRDIRRN DRGPGGPRGI VGGGMMRDRD GRGPPRGGM AQKLGSGRGT

image



Background

Alternative Names Ras GTPase-activating protein-binding protein 2 (G3BP-2) (GAP SH3 domain-binding protein 2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human G3BP2 Protein
Catalog #: PEH5429

