


Product Name: Recombinant Human CDK11A Protein
Catalog #: PEH5435



Summary

Name	CD11A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UQ88
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTLDEILQEK KRRKEQEEKA EIKRLKNSDD RDSKRDSLEE GELRDHCMEI TIRNSPYRRE DSMEDRGEED DSLAIKPPQQ MSRKEKVHHR KDEKRKEKWK HARVKEREHE

image



Background

Alternative Names	Cyclin-dependent kinase 11A (EC 2.7.11.22) (Cell division cycle 2-like protein kinase 2) (Cell division protein kinase 11A) (Galactosyltransferase-associated protein kinase p58/GTA) (PITSLRE serine/threonine-protein kinase CDC2L2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CDK11A Protein
Catalog #: PEH5435

