


Product Name: Recombinant Human ZBTB1 Protein
Catalog #: PEH5438



Summary

Name	ZBTB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2K1
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CGFGFSCEKL LDEHVLCTN RHLYQNTRSY HRIVDIRDGK DSNIAEFGE KDSSKTFSQAQ TDKYRGDTSQ AADDASTTG SRKSSTVESE IASEEKSRAA ERKRRIIKME PEDIPTDELK DFNIIKVTDK DCNESTDNDE LEDEPEEPFY RYYVEEDVSI KKSGRKTLP RMSVSADERG GLENMRPPNN SSPVQEDAEN ASCELCGLTI TEEDLSSHLY AKHIENICAC GKCGQILVKG RQLQEHAQRC GEPQDLTMNG LGNTEEKMDL EENPDEQSEI RDMFVEMLDD FRDNHYQINS IQKKQLFKHS ACPFRCPNCG QRFETENLVV EHMSSCLDQD MFKSAIMEEN

image



Background

Alternative Names Zinc finger and BTB domain-containing protein 1

Product Name: Recombinant Human ZBTB1 Protein
Catalog #: PEH5438



Note

For Research Use Only , Not for Diagnostic Use