


Product Name: Recombinant Human ZKSCAN5 Protein
Catalog #: PEH5439



Summary

Name	ZKSC5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2L8
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLDPPEAETSQ EQEDLFIVKV EEEDCTWMQE YNPPTFETFY QFRHFQYHE ASGPREALSQ LRVLCCEWLR PELHTKEQIL ELLVLEQFLT ILPEEFQPWV REHHPESGEE AVAVIENIQR ELEERRQQIV ACPDVLPRKM ATPGAVQESC SPHPLTVDTQ PEQAPQKPRL LEENALPVLQ VPSLPLKDSQ ELTASLLSTG SQKLVKIEEV ADVAVSFILE EWGHLDQSQK SLYRDDRKEN YGSITSMGYE SRDNMELIVK QISDDSESHW VAPEHTERS SV PQDPDFAEVS DLKGMVQRWQ VNPTVGKSRQ NPSQKRDLD A

image



Background

Alternative Names	Zinc finger protein with KRAB and SCAN domains 5 (Zinc finger protein 95 homolog) (Zfp-95)
--------------------------	--

Product Name: Recombinant Human ZKSCAN5 Protein
Catalog #: PEH5439



Note

For Research Use Only , Not for Diagnostic Use