


Product Name: Recombinant Human KLF13 Protein
Catalog #: PEH5441



Summary

Name	KLF13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2Y9
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGAAAAPPSP AWSEPEPEAG LEPEREPGPA GSGEPGLRQR VRRGRSRADL ESPQRKHKCH YAGCEKVYGK SSKLKAHLRT HTGERPFACS WQDCNKKFAR SDELARHYRT

image



Background

Alternative Names	Krueppel-like factor 13 (Basic transcription element-binding protein 3) (BTE-binding protein 3) (Novel Sp1-like zinc finger transcription factor 1) (RANTES factor of late activated T-lymphocytes 1) (RFLAT-1) (Transcription factor BTEB3) (Transcription factor NSLP1)
--------------------------	---

Note

Product Name: Recombinant Human KLF13 Protein
Catalog #: PEH5441



For Research Use Only , Not for Diagnostic Use