


Product Name: Recombinant Human HDGFL3 Protein
Catalog #: PEH5445



Summary

Name	HDGR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y3E1
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KFGKSNKRKG FNEGLWEIEN NPGVKFTGYQ AIQQSSSET EGEGGNTADA SSEEEGDRVE EDGKGKRKNE KAGSKRKKSYS TSKKSSKQSR KSPGDEDDKD CKEEENKSSS EGGDAGNDTR NTTSDLQKTS

image



Background

Alternative Names	Hepatoma-derived growth factor-related protein 3 (HRP-3) (Hepatoma-derived growth factor 2) (HDGF-2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use