


Product Name: Recombinant Human HECTD4 Protein
Catalog #: PEH5451



Summary

Name	HECD4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y4D8
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SICSSQGISQ TVSDLSVDPL PAGLELPIPP GLLEPHAVSS QESLDISLCS TGSLGSLGSL GEPLDNAETA SVSDMGSMYV VTSLDNQPLA ARPIKGFVV EIRSRAKIEK IRASLFNNND LIGLSSLDGE DELMEMSTEE

image



Background

Alternative Names	Probable E3 ubiquitin-protein ligase HECTD4 (EC 2.3.2.26) (HECT domain-containing protein 4) (HECT-type E3 ubiquitin transferase HECTD4)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use