


Product Name: Recombinant Human AFG3L2 Protein
Catalog #: PEH5452



Summary

Name	AFG32_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y4W6
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVLQPEEKKT VAYHEAGHAV AGWYLEHADP LLKVSIIPRG KGLGYAQYLP KEQYLYTKEQ LLDRMCMTLG GRVSEEIFFG RITTGAQDDL RKVTSAYAQ IVQFGMNEKV GQISFDLPRQ GDMVLEKPYS EATARLIDDE VRILINDAYK RTVALLTEKK ADVEKVALLL LEKEVLDKND MVELLGPRPF AEKSTYEETFV EGTGSLDEDT SLPEGLKDOWN KEREKEKEEP

image



Background

Alternative Names AFG3-like protein 2 (EC 3.4.24.-) (Paraplegin-like protein)

Note

Product Name: Recombinant Human AFG3L2 Protein
Catalog #: PEH5452



For Research Use Only , Not for Diagnostic Use