


Product Name: Recombinant Human PAXBP1 Protein
Catalog #: PEH5453



Summary

Name	PAXB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5B6
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELNSSAESEQ PLDKTGHVKD TNQEDGVIIS EHGEDEM DME SEKEEEKPKT GGAFSNALSS LNVLRPGEIP DAAFIHAARK KRQMARELGD FTPHDNEPGK GRLVREDEND ASDDEDDDEK RRIVFSVKEK SQRQKIAEEI GIEGSDDDAL

image



Background

Alternative Names PAX3- and PAX7-binding protein 1 (GC-rich sequence DNA-binding factor 1)

Note

For Research Use Only , Not for Diagnostic Use