


Product Name: Recombinant Human PCDHGB7 Protein
Catalog #: PEH5456



Summary

Name	PCDGJ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5F8
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LREDVPPGTS ILRVKATDQD EGINSEITYS FFGVADKAQH VFSLDYTTGN ILTQQPLDFE EVERYTINIE AKDRGSLSTR CKVIVEVVDE NDNSPEIIIT SLSDQIMEDS PPGVVVALFK TRDQDSGENG EVRCSLSRGV PFKIHSSSNN YYKLVTDEAL DREQTPEYNV TIAATDRGKP PLSSSKTITL HITDVNDNAP VFGQSAYLVH VPENNQPGAS IAQVSASDPD FGLNGRVSY S LIASDLESRT LSSYVSVAQ SGVVFAQRAF DHEQLRTFEL TLQARDQGSP ALSANVSLRV LVGDRNDNAP

image



Background

Alternative Names Protocadherin gamma-B7 (PCDH-gamma-B7)

Product Name: Recombinant Human PCDHGB7 Protein
Catalog #: PEH5456



Note

For Research Use Only , Not for Diagnostic Use