


Summary

Name	PCDGI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5F9
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSENLPPGSP VLQVTATDQD EGVNAEINYY FRSTAQSTKH MFSLDEKTGM IKNNQSDFE DVERYTMEVE AKDGGGLSTQ CKVIIILDE NDNSPEIIIT SLSDQILENS PPGMVVALFK TRDLDFGGNG EVRCNIETDI PFKIYSSSNN YYKLVDGAL DREQTPEYNV TIVATDRGKP PLSSRSITL YVADINDNAP VFDQTSYVVH VAENPPGAS IAQVSASDPD LGLNGHISYS IVASDLEPLA VSSYVSVSAQ SGVVFAQRAF DHEQLRAFAL TLQARDHGSP TLSANVSLRV LVGDRNDNAP

image



Background

Alternative Names Protocadherin gamma-B6 (PCDH-gamma-B6)

Product Name: Recombinant Human PCDHGB6 Protein
Catalog #: PEH5457



Note

For Research Use Only , Not for Diagnostic Use