


Summary

Name	PCDG3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5H0
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VWENVPVGTR LLTVNATDPD EGFNAQVSYI LDKMPGKIAE IFHLNSVSGE VSILKSLDYE DAMFYEIKIE AQDGPGLLSR AKILVTVLDV NDNAPEITIT SLTSSVPEEG TVGREIALID VHDRDSGQNG QVEVFVLGNL PFKLEKSIDQ YYRLVTATSL DREQISEYNI SLRASDGGSP PLSTETHITL HVIDINDNPP TFPHLSYSAY IPENNPRGAS IFSVTAQDPD SNNNARITYA LTEDTLQGAP LSSFVSINSN TGVLVALRSF DYEQFRDLKL LVTASDSGNP PLSSNVSLNL FVLDQNDNAP

image



Background

Alternative Names Protocadherin gamma-A3 (PCDH-gamma-A3)

Product Name: Recombinant Human PCDHGA3 Protein
Catalog #: PEH5459



Note

For Research Use Only , Not for Diagnostic Use