


Product Name: Recombinant Human PCDHA3 Protein
Catalog #: PEH5463



Summary

Name	PCDA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5H8
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLENAPNGTL VVTVNATDLD EGVNKDIAYS FNTDMSADIL SKFHLDPVNG QISVKGNIDF EESKSYEIQV EATDKGNPPM SDHCTVLLEI VDINDNVPEL VIQSLSLPVL EDSPLSTVIA LISVSDRDSG VNGQVTCSLT PHVPFKLVST FKNYYSLVLD SPLDRESVSA YELVVTARDG GSPSLWATAS VSVEVADVND NAPAFSQSEY TVFVKENNPP GCHIFTVSAR DADAQENALV SYSLVERRVG ERALSSYVSV HAESGKYYAL QPLDHEEEL LQFQVSARDA GVPPLGSNVT LQVFVLDEND

image



Background

Alternative Names	Protocadherin alpha-3 (PCDH-alpha-3)
--------------------------	--------------------------------------

Product Name: Recombinant Human PCDHA3 Protein
Catalog #: PEH5463



Note

For Research Use Only , Not for Diagnostic Use