


Product Name: Recombinant Human PCDHAC2 Protein
Catalog #: PEH5464



Summary

Name	PCDC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5I4
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YSPVPENATP NTIVAVLSVN DQDSGPNRKV SLGLEATLPF RLNGFGNSYT LVVSGPLDRE RVAVYNITVT ATDGGIPQLT SLRTLKVEIS DINDNPPSFL EDSYSIIQEQ NNLPGVLLCT VQATDPDEKE NAEVTYSLLE REIQGLPVTS YVSINSASGS LYAVNSFDYE KFREFVTVTE AQDKGSPPLS STVTANVYVV DMNDHAPHIL YPTSTNSSAA

image



Background

Alternative Names Protocadherin alpha-C2 (PCDH-alpha-C2)

Note

Product Name: Recombinant Human PCDHAC2 Protein
Catalog #: PEH5464



For Research Use Only , Not for Diagnostic Use