


Summary

Name	UTP18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5J1
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IAVAAEEER RLRQRNRLRL EEDKPAVERC LEELVFGDVE NDEDALLRRL RGRVQEHED SGDSEVENEA KGNFPPQKKP VVWDEEDEDE EMVDMMNRF RKDMMKNASE SKLSKDNLKK RLKEEFQHAM GGVPWAETT KRKTSSDES EEDEDLLQR TGNFISTSTS LPRGILKMKN CQHANAERPT VARISSVQFH PGAQIVMVAG LDNAVSLFQV DGKTNPKIQS IYLERFPIFK ACFSANGEEV LATSTHSKVL

image



Background

Alternative Names	U3 small nucleolar RNA-associated protein 18 homolog (WD repeat-containing protein 50)
--------------------------	--

Product Name: Recombinant Human UTP18 Protein
Catalog #: PEH5465



Note

For Research Use Only , Not for Diagnostic Use