


Product Name: Recombinant Human MAFB Protein
Catalog #: PEH5466



Summary

Name	MAFB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5Q3
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PHHAYPGAGV AHDELGPHAH PHHHHHHQAS PPPSSAASPA QQLPTSHPGP GPHATASATA AGGNGSVEDR FSDDQLVSMS VRELNRHLRG FTKDEVIRLK QKRRTLKNRG YAQSCRYKRV QQKHLENEK TQLIQQVEQL KQEVSLARE RDAYKVKCEK LANSGFREAG STDSPSSPE

image



Background

Alternative Names	Transcription factor MafB (Maf-B) (V-maf musculoaponeurotic fibrosarcoma oncogene homolog B)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MAFB Protein
Catalog #: PEH5466

