


Product Name: Recombinant Human ARFGEF2 Protein
Catalog #: PEH5468



Summary

Name	BIG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6D5
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQQKEIIIEHG IELFNKKPKR GIQFLQEQGM LGTSVEDIAQ FLHQEERLDS TQVGDFLGDS ARFNKEVMYA YVDQLDFCEK EFVSALRTFL EGFRLPGEAQ KIDRLMEKFA ARYIECNQGQ TLFASADTAY VLAYSIIMLT TDLHSPQVKN KMTKEQYIKM NRGINDSKDL PEEYLSSIYE EIEGKKIAMK

image



Background

Alternative Names	Brefeldin A-inhibited guanine nucleotide-exchange protein 2 (Brefeldin A-inhibited GEP 2) (ADP-ribosylation factor guanine nucleotide-exchange factor 2)
--------------------------	--

Note

Product Name: Recombinant Human ARFGEF2 Protein
Catalog #: PEH5468



For Research Use Only , Not for Diagnostic Use