


Product Name: Recombinant Human CDY2A Protein
Catalog #: PEH5469



Summary

Name	CDY2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6F7
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IVDKRQDKNG NTQYLVRWKG YDKQDDTWEP EQHLMNCEKC VHDFNRRQTE KQKKLTWTTT SRIFSNNARR RTSRSTKANY SKNSPKTPVT DKHHRSKNCK LFAASKNVRR KAASTLSDTK NMEIINSTIE TLAPDSPFDH KKTVSGFQKL EKLDPIAADQ QDTVVFVKVTE GKLLRDPLSH PGAEQTGIQN KTQMHPLMSQ MSGSVTASMA

image



Background

Alternative Names Testis-specific chromodomain protein Y 2 (EC 2.3.1.48)

Note

Product Name: Recombinant Human CDY2A Protein
Catalog #: PEH5469



For Research Use Only , Not for Diagnostic Use