



Product Name: Recombinant Human ZFP37 Protein
Catalog #: PEH5471



Summary

Name	ZFP37_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6Q3
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKILKQNLDL PDHSRNCVKR KSDAAKEHKK SFNHSLSDTR KGKKQTGKKH EKLSSHSSSD KCKNTGKKHD KLCCHSSSHI KQDKIQTGEK HEKSPSLSSS TKHEKPQACV KPYECNQCGK VLSHKQGLID HQRVHTGEKP YECNECGIAF SQKSHLVVHQ RTHTGEKPYE

image



Background

Alternative Names Zinc finger protein 37 homolog (Zfp-37)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZFP37 Protein
Catalog #: PEH5471

