


Product Name: Recombinant Human PGF Protein
Catalog #: PEH5476



Summary

Name	PLGF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P49763
Host	E.coli
Species	Human
Molecular Mass	65 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPAVPPQQWA LSAGNGSSEV EVVPFQEVWG RSYCRALERL VDVVSEYPSE VEHMFSPSCV SLLRCTGCCG DENLHCVPVE TANVTMQLLK IRSGDRPSYV ELTFSQHVRC ECRHSPGRQS PDMPGDFRAD APSFLPPRRS LPMLFRMEWG CALTGSQSAV WPSSPVPEI PRMHPGRNGK KQQRKPLREK MKPERCGDAV PRR

image



Background

Alternative Names Placenta growth factor (PIGF)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PGF Protein
Catalog #: PEH5476

