

Product Name: Recombinant Human OSMR Protein
Catalog #: PEH5485



Summary

Name	OSMR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q99650
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEEERIAGTE GGFSLSWKPQ PGDVIGYVVD WCDHTQDVLG DFQWKNVGP TTSTVISTDA FRPGVRYDFR IYGLSTKRIA CLEKKTGYS QELAPSDNPH VLVDTLTSHS FTLSWKDYST ESQPGFIQGY HVYLKSKARQ CHPRFEKAVL SDGSECCKYK IDNPEEKALI VDNLPESFY EFFITPFTSA GEGPSATFTK VTTPDEHSSM

image



Background

Alternative Names	Oncostatin-M-specific receptor subunit beta (Interleukin-31 receptor subunit beta) (IL-31 receptor subunit beta) (IL-31R subunit beta) (IL-31R-beta) (IL-31RB)
--------------------------	--

Product Name: Recombinant Human OSMR Protein
Catalog #: PEH5485



Note

For Research Use Only , Not for Diagnostic Use